

Results from Compound Optimization



Compound Optimization in MS and MS/MS

```

14:03:11: Optimizing S-Lens RF Amplitude for ion 340.18 m/z
14:03:17: Previous Setting = 149.70, New Setting = 86.99
14:03:17: Maximum Intensity = 9.14e+06
14:03:17: 56 % Improvement
14:03:18: Old Parent Mass: 340.176, New Parent Mass: 340.283
14:03:18: Optimizing collision energy at 1.5 mTorr
14:03:18: Waiting for the collision gas to stabilize
14:03:39: Finding the product ions of ion 340.3 m/z
14:04:17: Constructing the breakdown curve of ion 340.3 m/z
14:04:19: Product Ion: 251.16 Maximum Intensity: 2.18e+07
14:04:21: Product Ion: 162.20 Maximum Intensity: 3.48e+06
14:04:23: Product Ion: 179.06 Maximum Intensity: 4.70e+06
14:04:24: Product Ion: 193.07 Maximum Intensity: 4.43e+06
14:04:26: Product Ion: 104.12 Maximum Intensity: 3.98e+06
14:04:28: Product Ion: 180.07 Maximum Intensity: 2.08e+06
14:04:30: Product Ion: 158.11 Maximum Intensity: 1.80e+06
14:04:32: Product Ion: 209.11 Maximum Intensity: 1.25e+06
14:04:32: Collision Energy Optimization Results:
14:04:32: Product Ions (m/z) Coll. Energy (v)
14:04:32: 251.16 19
14:04:32: 162.20 20
14:04:32: 179.06 35
14:04:32: 193.07 33
14:04:32: 104.12 32
14:04:32: 180.07 40
14:04:32: 158.11 39
14:04:32: 209.11 37
14:04:33: Finish compound optimization
    
```

Comments:


Signature: \_\_\_\_\_




Signature: \_\_\_\_\_

# TSQ Tune - Tune and Calibration

SRIGAMP DAC Scan #A: 19 Peak 85.72; 9.21e+006



S#: 22541 FULL: PRO: 340 CT: 0.53 #A: 9 1.74e5



Signature: \_\_\_\_\_