

Results from Compound Optimization



Compound Optimization in MS and MS/MS

```

14:05:15: Optimizing S-Lens RF Amplitude for ion 340.18 m/z
14:05:21: Previous Setting = 149.70, New Setting = 92.68
14:05:21: Maximum Intensity = 1.02e+07
14:05:21: 29 % Improvement
14:05:22: Old Parent Mass: 340.176, New Parent Mass: 340.279
14:05:22: Optimizing collision energy at 1.5 mTorr
14:05:22: Waiting for the collision gas to stabilize
14:05:43: Finding the product ions of ion 340.3 m/z
14:06:21: Constructing the breakdown curve of ion 340.3 m/z
14:06:23: Product Ion: 251.17 Maximum Intensity: 2.35e+07
14:06:25: Product Ion: 179.06 Maximum Intensity: 5.46e+06
14:06:26: Product Ion: 162.19 Maximum Intensity: 6.18e+06
14:06:28: Product Ion: 193.08 Maximum Intensity: 3.42e+06
14:06:30: Product Ion: 104.13 Maximum Intensity: 4.07e+06
14:06:32: Product Ion: 180.07 Maximum Intensity: 2.04e+06
14:06:34: Product Ion: 158.11 Maximum Intensity: 1.80e+06
14:06:35: Product Ion: 144.10 Maximum Intensity: 1.63e+06
14:06:36: Collision Energy Optimization Results:
14:06:36: Product Ions (m/z) Coll. Energy (v)
14:06:36: 251.17 20
14:06:36: 179.06 36
14:06:36: 162.19 21
14:06:36: 193.08 36
14:06:36: 104.13 32
14:06:36: 180.07 35
14:06:36: 158.11 43
14:06:36: 144.10 48
14:06:37: Finish compound optimization
    
```

Comments:


Signature: _____




Signature: _____

TSQ Tune - Tune and Calibration

SRIGAMP DAC Scan #A: 19 Peak 92.68; 1.08e+007



S#: 22953 FULL: PRO: 340 CT: 0.52 #A: 9 2.44e5



Signature: _____