

Results from Compound Optimization



Compound Optimization in MS and MS/MS

```

14:00:59: Optimizing S-Lens RF Amplitude for ion 340.18 m/z
14:01:05: Previous Setting = 149.70, New Setting = 137.28
14:01:05: Maximum Intensity = 6.31e+06
14:01:05: 10 % Improvement
14:01:06: Old Parent Mass: 340.176, New Parent Mass: 340.286
14:01:06: Optimizing collision energy at 1.5 mTorr
14:01:06: Waiting for the collision gas to stabilize
14:01:27: Finding the product ions of ion 340.3 m/z
14:02:05: Constructing the breakdown curve of ion 340.3 m/z
14:02:07: Product Ion: 251.17 Maximum Intensity: 2.15e+07
14:02:09: Product Ion: 179.07 Maximum Intensity: 4.11e+06
14:02:10: Product Ion: 193.08 Maximum Intensity: 3.89e+06
14:02:12: Product Ion: 104.14 Maximum Intensity: 3.64e+06
14:02:14: Product Ion: 162.20 Maximum Intensity: 4.18e+06
14:02:16: Product Ion: 158.10 Maximum Intensity: 1.52e+06
14:02:18: Product Ion: 144.08 Maximum Intensity: 1.48e+06
14:02:19: Product Ion: 73.14 Maximum Intensity: 1.16e+06
14:02:20: Collision Energy Optimization Results:
14:02:20: Product Ions (m/z) Coll. Energy (v)
14:02:20: 251.17 19
14:02:20: 179.07 34
14:02:20: 193.08 30
14:02:20: 104.14 32
14:02:20: 162.20 20
14:02:20: 158.10 40
14:02:20: 144.08 53
14:02:20: 73.14 31
14:02:21: Finish compound optimization
    
```

Comments:


Signature: _____




Signature: _____

TSQ Tune - Tune and Calibration

SRIGAMP DAC Scan #A: 19 Peak 135.70; 6.56e+006



S#: 22119 FULL: PRO: 340 CT: 0.51 #A: 9 2.02e5



Signature: _____