

Results from Compound Optimization



Compound Optimization in MS and MS/MS

```

13:40:15: Optimizing S-Lens RF Amplitude for ion 336.18 m/z
13:40:20: Previous Setting = 102.39, New Setting = 108.18
13:40:20: Maximum Intensity = 3.74e+06
13:40:20: 2 % Improvement
13:40:21: Old Parent Mass: 336.176, New Parent Mass: 336.272
13:40:21: Optimizing collision energy at 1.5 mTorr
13:40:21: Waiting for the collision gas to stabilize
13:40:42: Finding the product ions of ion 336.3 m/z
13:41:21: Constructing the breakdown curve of ion 336.3 m/z
13:41:22: Product Ion: 247.15 Maximum Intensity: 4.06e+06
13:41:24: Product Ion: 179.07 Maximum Intensity: 8.19e+05
13:41:26: Product Ion: 191.06 Maximum Intensity: 9.37e+05
13:41:28: Product Ion: 158.19 Maximum Intensity: 8.79e+05
13:41:30: Product Ion: 102.11 Maximum Intensity: 6.09e+05
13:41:31: Product Ion: 156.09 Maximum Intensity: 3.70e+05
13:41:33: Product Ion: 205.08 Maximum Intensity: 4.53e+05
13:41:35: Product Ion: 144.09 Maximum Intensity: 4.53e+05
13:41:35: Collision Energy Optimization Results:
13:41:35: Product Ions (m/z) Coll. Energy (v)
13:41:35: 247.15 18
13:41:35: 179.07 33
13:41:35: 191.06 32
13:41:35: 158.19 22
13:41:35: 102.11 32
13:41:35: 156.09 42
13:41:35: 205.08 34
13:41:35: 144.09 52
13:41:37: Finish compound optimization
    
```

Comments:


Signature: _____




Signature: _____

TSQ Tune - Tune and Calibration

SRIGAMP DAC Scan #A: 19 Peak 107.86; 3.63e+006



S#: 19490 FULL: PRO: 336 CT: 0.61 #A: 9 9.67e4



Signature: _____